



DENTAL IMPLANTS

MADE IN GERMANY



PRODUCT CATALOG

SINCE 1995



ta-implants.com



TA Switch / Plus
The Innovative Implant System




PREMIUM QUALITY

About us

Founded in 1995, TA Dental Implants embodies Germany's renowned engineering tradition and precision manufacturing in the field of dental implantology. Since our inception, we have been committed to delivering innovative solutions, exceptional production quality, and premium-level reliability to dental professionals and patients worldwide.


Our entire production process takes place in state-of-the-art facilities in Germany, where each TA Dental Implants product is meticulously crafted with advanced engineering and subjected to rigorous quality control. Designed to meet the highest standards of biocompatibility, durability, and aesthetics, our implants reflect our unwavering dedication to excellence.

At TA Dental Implants, we continuously invest in scientific research and cutting-edge technology to develop pioneering solutions that set new benchmarks in the industry. With our unwavering commitment to superior German quality, we take pride in being a trusted partner in creating healthy, confident smiles.



Content

Dental Implants	3-6
TA Switch / Plus Passive	3
TA Switch / Plus Aggressive	5
Packaging	7-8
Surgery	9-14
Surgical Drills & Taps	9
Drilling Protocol	11
Surgical Instruments	13
Prosthetic Parts	15-30
Healing Abutments & Cover Screws	15
Impression Components	17
Prosthetic Profiles	19
Prosthetic ClicLoc	23
All-On-TA System	25
CAD-CAM Parts	29



TA Switch / Plus Passive

22 ° Conical Connection

Color Code

- Ø 3.3
- Ø 3.8
- Ø 4.1
- Ø 4.8
- Ø 5.5

Micro Threads

Internal Hexagonal Indexing

2.1mm

Regular Platform


2.5mm

Wide Platform

Sinus Protective Apex

Cylindrical Design

Connective
Contour



Sandblasted and
Acid-Etched Surface

10000x

Surface roughness (Sa): 1.5 µm +/- 0.5 µm

implants



Ø 3.3

D: 3.3 mm - AD: 3.1 mm

Length	Product ID
8	2013308A
10	2013310A
12	2013312A
14	2013314A

Regular Platform



Ø 3.8

D: 3.8 mm - AD: 3.6 mm

Length	Product ID
8	2013808A
10	2013810A
12	2013812A
14	2013814A

Regular Platform



Ø 4.1

D: 4.1 mm - AD: 3.9 mm

Length	Product ID
8	2014108A
10	2014110A
12	2014112A
14	2014114A

Regular Platform



Ø 4.8

D: 4.8 mm - AD: 4.6 mm

Length	Product ID
8	2014808A
10	2014810A
12	2014812A
14	2014814A

Wide Platform




Ø 5.5

D: 5.5 mm - AD: 5.3 mm

Length	Product ID
8	2015508A
10	2015510A
12	2015512A
14	2015514A

Wide Platform



TA Switch / Plus **Passive - Short**

L - 6 mm



Ø 3.8



Ø 4.1



Ø 4.8



Ø 5.5

2013806A | 2014106A | 2014806A | 2015506A | Product ID

Regular Platform

TA Switch / Plus Aggressive

22 ° Conical Connection

Color Code

- Ø 3.3
- Ø 3.8
- Ø 4.1
- Ø 4.8
- Ø 5.5

Micro Threads

Internal Hexagonal Indexing

2.1mm Regular Platform


2.5mm Wide Platform

Deep Apical Threads

Sinus Protective Apex

Conical Design

Connective
Contour



Sandblasted and
Acid-Etched Surface

10000x

Surface roughness (Sa): 1.5 µm +/- 0.5 µm

implants




Length	Product ID
8	2023308A
10	2023310A
12	2023312A
14	2023314A

Regular Platform

Ø 3.3

D: 3.3 mm - AD: 2.7 mm



Length	Product ID
8	2023808A
10	2023810A
12	2023812A
14	2023814A

Regular Platform

Ø 3.8

D: 3.8 mm - AD: 3.2 mm



Length	Product ID
8	2024108A
10	2024110A
12	2024112A
14	2024114A

Regular Platform

Ø 4.1

D: 4.1 mm - AD: 3.5 mm



Length	Product ID
8	2024808A
10	2024810A
12	2024812A
14	2024814A

Wide Platform

Ø 4.8

D: 4.8 mm - AD: 4.2 mm




Length	Product ID
8	2025508A
10	2025510A
12	2025512A
14	2025514A

Wide Platform

Ø 5.5

D: 5.5 mm - AD: 4.9 mm



TA Switch / Plus Packaging

Aggressive
Implant Box



Passive
Implant Box




Prosthetic
Component
Package



Implant Blister



Remove the Tyvek film from the blister pack to reveal the implant container.
(Caution: Removing the film breaks the sterile barrier.)




How to Use The Implant Container




Implant Container




Slide the container cover to open and expose the implant.




Apply gentle pressure to the container to keep the implant stable while positioning the implant driver.



Ensure the implant driver is properly seated in the implant before releasing the pressure on the container, then remove the implant from the container.




The implant cover is located at the bottom of the implant container.



Remove the implant cover from the container by turning it counterclockwise with a hex driver.

Surgical Drills & Taps



Twist Drills

Step Drills

DLC Coating

First Step

Second Step



Surgical Drills & Taps



204001

Triangular Drill



204002

Ø 2.0 Pilot Drill



204047

Round Drill



204046

Bone Profiler Drill



205009

Drill Extension



204004

Ø 3.3 Twist Drill



204005

Ø 3.8 Twist Drill



204006

Ø 4.1 Twist Drill



204007

Ø 4.8 Twist Drill



204008

Ø 5.5 Twist Drill



204036

Ø 3.3 Step Drill



204037

Ø 3.8 Step Drill



204038

Ø 4.1 Step Drill



204039

Ø 4.8 Step Drill



204040

Ø 5.5 Step Drill



204009

Ø 3.3 Cortical Drill



204010

Ø 3.8 Cortical Drill



204011

Ø 4.1 Cortical Drill



204012

Ø 4.8 Cortical Drill



204013

Ø 5.5 Cortical Drill



204014

Ø 3.3 Tap Drill



204015

Ø 3.8 Tap Drill



204016

Ø 4.1 Tap Drill



204017

Ø 4.8 Tap Drill



204018

Ø 5.5 Tap Drill

Ø 3.3

Ø 3.8

Ø 4.1


Ø 4.8

Ø 5.5

Drilling Protocol

TA Switch / Plus Passive

Use Twist drills
specifically designed
for Passive implants.




Drilling Protocol

TA Switch / Plus

Aggressive

Use step drills
specifically designed
for Aggressive implants.



Surgical Instruments



Bone Profiler Drill removes excess bone around the implant site for all abutment types, ensuring accurate and stable seating of prosthetic components.



The ISO Key is an adapter that enables handpieces to be used with a torque ratchet and also serves as a torque wrench for straight All-on-TA abutments.



205001

Short Screwdriver



205002

Long Screwdriver



205003

Short Screwdriver HP



205004

Long Screwdriver HP



205005

Short Ball Screwdriver



205006

Long Ball Screwdriver



205007

Short Ball Screwdriver HP



205008

Long Ball Screwdriver HP



205012

Short Implant Driver




205013

Long Implant Driver



205016

Short Implant Driver HP



205017

Long Implant Driver HP

Regular Platform



205014

Short Implant Driver



205015

Long Implant Driver



205018

Short Implant Driver HP



205019

Long Implant Driver HP

Wide Platform



205020

Directional Pin



205023

20° Angled Pin



205024

30° Angled Pin



205025

Guide Pin Regular



205026

Guide Pin Wide



205028

Surgical Tray.



205011

Handle ISO Connection



205023

Dept Gauge



205010

Torque Ratchet

Healing Abutments

Color Code

● Regular Platform

● Wide Platform




Available in a wide range of diameters and lengths to suit all indications

Biological
S - Line
Design

1 mm

Perfect abutment fit



Prosthetic

Healing Abutments



D - Ø4.0

Gingival Height	Product ID
2 mm	203015
4 mm	203016
6 mm	203017

Regular Platform



D - Ø4.5

Gingival Height	Product ID
2 mm	203018
4 mm	203019
6 mm	203020

Regular Platform



D - Ø5.5

Gingival Height	Product ID
2 mm	203021
4 mm	203022
6 mm	203023

Regular Platform



D - Ø5.5

Gingival Height	Product ID
2 mm	203024
4 mm	203025
6 mm	203026


Wide Platform



D - Ø6.5

Gingival Height	Product ID
2 mm	203027
4 mm	203028
6 mm	203029

Wide Platform



Cover Screws



203001

Regular Platform



203002

Wide Platform

Impression


Components



Color Code

● Regular Platform

● Wide Platform



Short and long impression posts to suit various clinical requirements.

It is compatible with the white closed tray transfer cap, regular platform closed tray impression post, and All-on-TA closed tray impression post.




It is compatible with the black closed tray transfer cap and wide platform closed tray impression posts.




Analogues are compatible with both conventional and digital impressions.

Impression Components



206001A

Impression coping
open tray
Short / Regular
Incl. Screw.



206002A

Impression coping
open tray
Long / Regular
Incl. Screw.



206003

Screw for Impression coping
open tray
Short / Regular



206004

Screw for Impression coping
open tray
Long / Regular



206011A

Impression coping closed tray
Regular / Incl. Screw and Cover




206013

Cover for Impression coping
closed tray / Regular




206009

Implant analog
Regular




206005A

Impression coping
open tray
Short / Wide
Incl. Screw.



206006A

Impression coping
open tray
Long / Wide
Incl. Screw.



206007

Screw for Impression coping
open tray
Short / Wide



206008

Screw for Impression coping
open tray
Long / Wide



206012A

Impression coping closed tray
Wide / Incl. Screw



206014

Cover for Impression coping
closed tray / Wide



206010


Implant analog
Wide

Cement - Retained Abutments

Aesthetic Line


Color Code


- Regular Platform
- Wide Platform




Available in straight and 20° angled options, they offer optimal adaptability for diverse clinical indications — combining precision, stability, and aesthetic flexibility.

Designed to support and shape the natural gingival contour, ensuring an optimal emergence profile and superior soft tissue aesthetics.




Cement - Retained **Abutments** / Aesthetic Line

D
Ø5.0


207004A

Universal Abutment
Regular - Incl. Screw

D
Ø4.0


207006A

Aesthetic Abutment
Regular - Incl. Screw

D
Ø5.0

207007A

Aesthetic Abutment
Regular - Incl. Screw

D
Ø4.5


207008A

20° Edge Aesthetic Abutment
Regular - Incl. Screw

D
Ø4.5

207009A

20° Angled face Aesthetic Abutment
Regular - Incl. Screw

207044A

Neo-Hybrid Abutment
Regular - Incl. Screw

D
Ø3.5

207045A

Temporary Abutment
Regular - Incl. Screw

D
Ø3.5

207046A

Temporary Abutment
Non - Hex
Regular - Incl. Screw

207001

Abutment Screw
Regular


207003

Abutment Screw
Neo - Hybrid


Regular Platform


D
Ø6.0

207005A

Universal Abutment
Wide - Incl. Screw

D
Ø5.0

207010A

Aesthetic Abutment
Wide - Incl. Screw

D
Ø6.0

207011A

Aesthetic Abutment
Wide - Incl. Screw

D
Ø6.0

207012A

20° Edge Aesthetic Abutment
Wide - Incl. Screw


D
Ø6.0

207013A

20° Angled Face Aesthetic Abutment
Wide - Incl. Screw

D
Ø4.5

207047A

Temporary Abutment
Wide - Incl. Screw

D
Ø4.5

207048A

Temporary Abutment Non-Hex
Wide - Incl. Screw

207002

Abutment Screw
Wide

Wide Platform

Cement - Retained Abutments


Bio Line



Color Code


- Regular Platform
- Wide Platform

0 ° 20 °




Available in straight and 20° angled options, they offer optimal adaptability for diverse clinical indications — combining precision, stability, and aesthetic flexibility.

The biologically shaped S-Line design of the abutments ensures optimal soft tissue support and natural emergence profile — promoting long-term esthetics and peri-implant health.




esthetic
prosthesis

Cement - Retained **Abutments** / Bio Line



Gingival Height	Product ID
1.5 mm	207014A
3 mm	207015A
5 mm	207016A

Straight Abutment Regular
Incl. Screw




Gingival Height	Product ID
1.5 mm	207017A
3 mm	207018A
5 mm	207019A

Straight Abutment Regular
Incl. Screw




Gingival Height	Product ID
1.5 mm	207020A
3 mm	207021A
5 mm	207022A

Straight Abutment Regular
Incl. Screw




Gingival Height	Product ID
1.5 mm	207023A
3 mm	207024A
5 mm	207025A

20° Angled Abutment Regular
Incl. Screw



Gingival Height	Product ID
1.5 mm	207026A
3 mm	207027A
5 mm	207028A


20° Angled Abutment Regular
Incl. Screw



Gingival Height	Product ID
1.5 mm	207029A
3 mm	207030A
5 mm	207031A


20° Angled Abutment Regular
Incl. Screw

Regular Platform



Gingival Height	Product ID
1.5 mm	207032A
3 mm	207033A
5 mm	207034A

Straight Abutment Wide
Incl. Screw




Gingival Height	Product ID
1.5 mm	207035A
3 mm	207036A
5 mm	207037A

Straight Abutment Wide
Incl. Screw



Gingival Height	Product ID
1.5 mm	207038A
3 mm	207039A
5 mm	207040A

20° Angled Abutment Wide
Incl. Screw



Gingival Height	Product ID
1.5 mm	207041A
3 mm	207042A
5 mm	207043A


20° Angled Abutment Wide
Incl. Screw

Wide Platform


Prosthetic ClickLoc



Each component is perfectly matched for long-term, comfortable use.



STANDARD INSERTS
Provides dual retention for strong stability and allows pivoting up to 20° of total divergence between two implants.



EXTENDED RANGE INSERTS
Allows pivoting movement to accommodate up to 40° of total divergence between two implants.



Prosthetic

Prosthetic ClicLoc



Gingival Height	Product ID
2 mm	208001
3.5 mm	208002
5 mm	208003

Regular Platform



Gingival Height	Product ID
2 mm	208004
3.5 mm	208005
5 mm	208006

Wide Platform



208008

Attachment Matrix black
reinforced retention:
10 – 13 N



208009

Attachment Matrix clear
divergence compensation:
< 20°
reinforced retention:
20.5 – 23.5 N



208010

Attachment Matrix pink
divergence compensation:
< 20°
reinforced retention:
12.5 – 15.5 N



208011

Attachment Matrix blue)
divergence compensation:
< 20°
reinforced retention:
7.5 – 10.5 N



208012

Attachment Matrix Green
Divergence Compensation:
20° – 40°
Reinforced Retention:
23.5 – 27.5 N



208013

Attachment Matrix Orange
Divergence Compensation:
20° – 40°
Reinforced Retention:
13 – 16.5 N



208014

Attachment Matrix Red
Divergence Compensation:
20° – 40°
Reinforced Retention:
5 – 8 N



208015

Titanium-Sleeve
incl. Matrix black



208016

Block-out Ring



208017

Impression post for
Clicloc System
incl. Matrix black



208018

Implant analog for
Clicloc System



208007

Locatorcap Al-On-TA



208019

Matrix Set-1
incl. Block-out-Ring
Divergence compensation:
< 20°



208020

Matrix Set-2
incl. Block-out-Ring
Divergence compensation:
20° – 40°



208021

Universal Tool
for Clicloc System



208022

Tool Clicloc for HP

Prosthetic All-on-TA




Color Code

- Regular Platform
- Wide Platform



Ergonomic carrier components allow for easy intraoral placement and precise positioning of the abutments.

All on TA abutments provide a standardized prosthetic platform with a 15° taper on one side. When combined with 30° angled abutments, they compensate for implant divergences of up to 90°.




Prosthetic




Gingival Height	Product ID
1.5 mm	209003
3 mm	209004
5 mm	209005

0° Abutment All on TA



GH(a)/GH(b)	Product ID
1.5 mm / 3.1 mm	209006A
3 mm / 4.7 mm	209007A

20° Angled Abutment All on TA Incl. Screw



GH(a)/GH(b)	Product ID
1.5 mm / 3.9 mm	209008A
3 mm / 5.4 mm	209009A

30° Angled Abutment All on TA Incl. Screw



209001

Angled Abutment Screw All-on-TA
Regular

Regular Platform



Gingival Height	Product ID
1.5 mm	209010
3 mm	209011
5 mm	209012

0° Abutment All on TA



GH(a)/GH(b)	Product ID
1.5 mm / 3.1 mm	209013A
3 mm / 4.7 mm	209014A

20° Angled Abutment All on TA Incl. Screw



GH(a)/GH(b)	Product ID
1.5 mm / 3.9 mm	209015A
3 mm / 5.4 mm	209016A



30° Angled Abutment All on TA Incl. Screw



209002

Angled Abutment Screw All-on-TA
Regular

Wide Platform



Prosthetic All-on-TA




Prosthetic



All on TA offers digital, open, and closed impression options — for ultimate flexibility and a smooth prosthetic workflow.

Perfect fit and precision-engineered components eliminate screw loosening and fracture — ensuring long-term stability and confidence.





209017

Okklusal Screw All-on-TA



209018A

Impression Coping
Open Tray for All-on-TA
Incl. screw



209019

Screw for Impression Coping
Open Tray All-on-TA



209020A

Impression Coping
Closed Tray for All-on-TA
Incl. screw and Cover



206013

Cover for Impression Coping
Closed Tray



209022

Implant Analog
All-on-TA



209023

Gingiva Sleeve Short
All-on-TA



209024

Gingiva Sleeve Long
All-on-TA



206025A

Gingiva Sleeve Long
All-on-TA / Peek



209026A

Scanadapter
All-on-TA
Incl. screw



209027A

Bridgecap Titanium
All-on-TA
Incl. screw



209028A

Bridgecap Castable
All-on-TA
Incl. screw



209029

Technical Screw
All-on-TA



209030

Insertion Instrument
All-on-TA
for Angled Abutment



209032

Temporary Cap
All-on-TA



209031

Drillbar All-on-TA

Prosthetic Cad-Cam



Both intraoral and extraoral scanning options are available. Intraoral scanning components are textured and color-marked to ensure greater accuracy and precision during scanning.

TA Flex abutments offers a reliable screw-retained restoration solution by allowing angulated screw access of up to 25°, effectively relocating the screw channel away from the esthetic zone.



Max 25°

Prosthetic

Prosthetic Cad-Cam




210003A

CAD-CAM Bonding Base
Plateau Ø 3.9 mm
Incl. Screw



210004A

CAD-CAM Bonding Base
Plateau Ø 4.5 mm
Incl. Screw



210009A

TA-Flex Regular



210006A

CAD-CAM Bonding Base
Non-Hex Plateau Ø 3.9 mm
Incl. Screw



210007A

CAD-CAM Bonding Base
Non-Hex Plateau Ø 4.5 mm
Incl. Screw



210011A

TA-Flex Non-Hex Regular



210015A

Intraoral Scan Body
Regular



210017A

Laboratory Scan Body
Regular



210013A


Pre-mil Abutment
Regular



210001


Ball-Screw TA-Flex Wide

Regular Platform




210005A

CAD-CAM Bonding Base
Plateau Ø 5.1 mm
Incl. Screw




210008A

CAD-CAM Bonding Base
Non-Hex Plateau Ø 5.1 mm
Incl. Screw



210010A

TA-Flex Wide



210012A

TA-Flex Non-Hex Wide




210016A

Intraoral Scan Body
Wide



210018A

Laboratory Scan Body
Wide



210014A

Pre-mil Abutment
Wide



210002

Ball-Screw TA-Flex Wide

Wide Platform



The Innovative Implant System



PREMIUM QUALITY

TA Zahnimplantate GmbH

Beim Hochwachturm 15
71332 Waiblingen
GERMANY

E-Mail : info@ta-implants.com